pro-fixTM
PRECISION FIXATION SYSTEM

INSTRUCTIONS FOR USE

EN
English

DE
Deutsch

FR
Français

ES
Español

IT
Italiano

PT
Português

SV
Svenska

DA
Dansk

NL
Nederlands

Introduction:

The user of Osteogenics Biomedical products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Osteogenics Biomedical disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use or installation of Osteogenics Biomedical products.

IMPORTANT: Read this entire package insert prior to use and follow all instructions carefully.

Description:

The *Pro-Fix*TM Precision Fixation System consists of screws of various sizes, the instrumentation that is used to place the screws, and the sterilization container.

Indications:

The *Pro-Fix*TM Precision Fixation System is used to stabilize, fixate, and/or support bone grafts, bone filling materials, and/or barrier membranes used for regeneration of bone in the oral cavity.

Instructions for Use:

A variety of surgical techniques may be used during placement and fixation of the devices. The placement sites and techniques used are at the discretion of the surgeon. Ensure that the appropriate instrumentation is available prior to surgery. Instrumentation includes *Pro-Fix*TM blades that fit the *Pro-Fix*TM handle or contra-angle hand pieces. These blades specifically fit all of the *Pro-Fix*TM screws.

Note: Use *Pro-Fix*TM screwdriver only with *Pro-Fix*TM screws.

Note: To preserve the surface finish, these implants should be handled with clean instruments dedicated to titanium, or a talc-free gloved hand.

Note: To disengage the *Pro-Fix*TM screwdriver from the screw tilt the screwdriver at a 45 degree angle to loosen and release the friction-fit connection between the driver and screw.

Note: when using *Pro-Fix*TM blades designed to fit contra-angle hand pieces, use the manufacturer's speed recommendations to produce the best results.

CAUTION: THIS PRODUCT IS NOT INTENDED TO BE THE SOLE MEANS OF SUPPORT. NO SUCH IMPLANT CAN WITHSTAND BODY LOADS WITHOUT THE SUPPORT OF BONE. IN THIS EVENT, BENDING, LOOSENING, DISASSEMBLY, AND/OR BREAKAGE OF THE IMPLANTS WILL EVENTUALLY OCCUR.

Material Specifications:

The screws are made of medical grade titanium alloy (ASTM F136), and are expressly warranted as being fabricated from the foregoing material specifications.

Cleaning and Sterilization Instructions:

The fixation screws are provided clean and **NON-STERILE** (requiring sterilization), and are able to withstand multiple sterilizations. The fixation screws are for single patient use only. The instrumentation (screwdriver blade and handle) is also provided clean and **NON-STERILE** (requiring sterilization), and is able to withstand multiple sterilizations. The instrumentation can be utilized repeatedly to insert the screws, but only after being cleaned and sterilized before every use. Manual cleaning is preferred. To clean, rinse blade and handle with warm running water, wash with appropriate detergent, and gently scrub with a nonmetallic brush. Both the instruments and screws must be sterilized in the *Pro-Fix*TM sterilization container prior to use. **CAUTION:** Do not exceed double wrapping to ensure adequate sterilization by the following validated parameters. The *Pro-Fix*TM fixation screws and instrumentation should be sterilized with moist heat, using either of the following validated steam sterilization guidelines.

Method	Cycle	Temperature	Exposure Time
Steam	Pre-vacuum	270°F 132°C	30 minutes
Steam	Gravity	250°F 121°C	60 minutes

How Supplied:

Unless provided sterile and clearly labeled as such, *Pro-Fix*TM *Precision Fixation System* screws are supplied **NON-STERILE**. They are intended for **SINGLE USE ONLY**. **NON-STERILE** products must be sterilized prior to use. Sterilize according to the suggested protocol described above. Do not use if package has been previously opened or damaged.

Patient Education:

It is the physician's responsibility to educate the patient and/or their representative(s) regarding oral-maxillofacial surgery. This should include a description of associated complications and an explanation of potential alternative products and treatments.

MRI Safety Information:

Non-clinical testing has demonstrated that *Profix Precision Fixation System* Screws are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T.
- Maximum spatial gradient of 3,000 gauss/cm (30 T/m).
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the *Profix Precision Fixation System* Screws are expected to produce a maximum temperature rise of less than 2.3°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device

extends approximately 3 mm from the implant when imaged with a gradient echo pulse sequence and a 3.0 T MR system.

Contraindications:

Contraindications include, but are not limited to, the following:

- Active infections.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be done to rule out this possibility prior to implantation.
- Conditions which tend to limit the patient's ability and/or willingness to cooperate and follow instructions during the healing period.
- Any degenerative disease, the progress of which could adversely affect the placement of the implant.
- Insufficient quantity of bone to securely anchor the implant.
- Blood supply limitations, radiation therapy, and/or previous infections which may tend to retard healing and increase the possibility of infection and/or rejection of the implant.
- Inadequate coverage with healthy tissue.
- Procedures in which a non-sterile environment exists, i.e. open cavities, such as sinuses.

Possible Adverse Effects:

All of the complications associated with surgery are possible. Complications and possible adverse effects associated with implants may, in addition, include the following:

- Decrease density of bone and/or necrosis of bone due to stress shielding.
- Vascular changes.
- Allergic reaction or metal sensitivity to the fixation devices.
- Nerve damage due to surgical trauma.
- Breakage of the fixation device due to a nonunion or delayed union of bony tissue.
- Bending or fracture of the fixation devices.
- Migration or loosening of the fixation devices.
- Pain, discomfort, and/or abnormal sensation due to presence of the fixation devices.
- Superficial and/or deep infection.
- Growth restriction.
- Passive transmigration of the fixation devices.
- Tissue staining.

Warnings and Precautions:

The physician is responsible for describing and explaining the following warnings, precautions and complications to the patient and/or his representative(s) prior to proceeding with the surgical procedures:

- Lint, fingerprints, talc and other surface contaminants or residues from latex gloves can cause foreign body or allergic reactions.
- Care must be taken to ensure that particulate contaminants are not introduced onto components during implantation or handling. This could result in improper performance of the system.
- An implanted device should never be re-used.
- Leftover implants which have been contaminated with blood or body fluids should be discarded as medical waste.
- Until bone healing is complete, the fixation provided by the system should be considered temporary and may not withstand extraordinary unsupported stresses.
- Detailed instructions on the use, limitations, and possible

- complications of the system should be given to the patient.
- Any decision to remove the system should take into consideration the potential risk to the patient of a second surgical procedure.
- Preoperative and operation procedures, including knowledge of surgical techniques and proper selection and placement of the implants, are important considerations in the successful utilization of the system.
- Radiation therapy has been shown to decrease the chance of successful result.
- Postoperatively as a precaution, before patients with implants receive any further surgery (such as dental procedures) prophylactic antibiotics may be considered, especially in high risk patients.


Warranty:

All products are warranted to be free of defects in material and workmanship. No warranty is made for any purpose other than in the product specifications and labeling.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician or dentist.

Labeling Symbols:

Symbols may be used on package labeling for easy identification.

-  Caution
-  Do Not Reuse
-  Manufacturer
-  Manufactured Date
-  Lot Number
-  Catalog Number
-  Not Sterile Product
-  Quantity
-  Authorized Representative in the European Community
-  Consult instructions for use
-  MR Conditional

 Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
www.osteogenics.com

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