CYTOPLAST[®]

RTMPLUG, RTMTAPE and RTMFOAM
Dental Wound Dressings

INSTRUCTIONS FOR USE

EN

English

PT

Português

DE

Deutsch

SV

Svenska

FR

Français

DA

Dansk

ES

Español

NL

Nederlands

IT

Italiano

EN

CYTOPLAST®

Dental Wound Dressings

INSTRUCTIONS FOR USE

Description:

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM are absorbent porous collagen matrices engineered from purified collagen derived from bovine dermis tissue. CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM are applied directly to the wound and protect the wound bed and delicate new tissue. CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM can be removed, replaced or left in situ. If left in situ the dressings will be essentially resorbed in 30 days. CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM are supplied sterile, non-pyrogenic and for single use only.

Indications and Usage:

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM are indicated for the management of moderate to severe oral wounds and sores, including:

- Denture sores
- Oral ulcers (non-infected or viral)
- Periodontal surgical wounds
- Suture sites
- Burns
- Extraction sites
- Surgical wounds
- Traumatic wounds

The wound dressings are designed for use in dental applications where, the use of the resorbable collagen wound dressing is considered beneficial by the clinician in managing the oral wound and the wound is not expected to heal on its own.

Contraindications:

This product is not indicated for patients with a known allergy to bovine-derived collagen materials. Large tissue defects should not be filled with several layers of Collagen Dental Wound Dressings.

Warning:

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM are not intended for use as a

bone void filler or bone grafting material.

Product Application:

The wound should be rinsed and excess fluid removed. CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM should be applied over the wound and held in place until adherence has been achieved. At the end of the procedure, CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM can be removed, replaced or left in situ.

Precautions:

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM should not be used on infected or contaminated wounds.

Discontinue the use of CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM and notify your doctor if excessive redness, pain, swelling or blistering occurs.

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM may be replaced if it falls off the wound prematurely. Replacement after premature removal of the Collagen Dental Wound Dressing is at the clinician's discretion.

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM have not been evaluated in pregnant women or children.

CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM cannot be re-sterilized or re-used. Open, unused CYTOPLAST® RTMPLUG, RTMTAPE and RTMFOAM must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

Clinicians should inform their patients prior to procedure that the dressing contains bovine derived collagen.

Bovine Dermis Collagen Dental Wound Dressings are not intended to be used on wounds that can heal spontaneously.

Adverse Reactions:

Hypersensitivity reactions have been noted with the use of other products containing collagen.

Storage:

The product should be stored at room temperature (15°C/59°F to 30°C/86°F). Avoid excessive heat and humidity.

How Supplied:

The dressings are supplied sterile and for single use only. They are available in a variety of shapes and sizes.

Individually packaged, ten (10) units per box:

Form	Size
RTMPLUG10	10 mm x 20 mm
RTMTAPE10	25 mm x 75 mm
RTMFOAM10	20 mm x 40 mm

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician or dentist.

Labeling Symbols:

Symbols may be used on some international package labeling for easy identification.

Manufacturer

Use by

Do Not Reuse

Caution

Lot Number

Catalog Number

Do not use if the product sterilization barrier or its packaging is compromised.

Temperature limitation

Sterilized Using Irradiation

Authorized Representative in the European Community

Federal (USA) law restricts this device to sale by or on the order of a physician or dentist

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Distributed by:

Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
www.osteogenics.com

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Collagen Matrix, Inc.
MDSS GmbH