Allotech, LLC – Cytoplast MicroDerm Package Insert

DESCRIPTION / USE: In accordance with FDA Article 21 CFR Part 1271 and Health Canada’s Safety of Human Cells, Tissues and Organs for Transplantation Regulations, this package contains donated human cells, tissues, and cellular and tissue-based products (HCT/Ps). Cytoplast™ MicroDerm is an acellular dermal matrix sourced from donated human tissue. Cytoplast™ MicroDerm is micro-surfaced, supplied in a range of sizes for surgical use by licensed clinicians, and is intended for homologous use only. Tissue is processed and packaged, solvent dehydrated, and terminally sterilized.

WARNINGS AND PRECAUTIONS: Human tissue has the potential to transmit infectious agents. Donor screening, processing treatments and laboratory testing follow stringent specifications to reduce the risk of infectious agent transmission. Do not use if the expiration date has been exceeded or if there is evidence of defects in package or label integrity. Do not use if the package integrity has been compromised. Once the user breaks the seal on the inner-most packaging, the tissue grafts must be transplanted or discarded. Do not sterilize or re-sterilize. It is the responsibility of the clinician to maintain tissue for transplantation according to recommended storage conditions. Do not use if tissue has not been stored according to the recommended STORAGE instructions. This tissue is intended for use by licensed clinicians.

Allografts are processed using some or all of the following agents: acetic acid, acetone, hydrogen peroxide, and sodium hydroxide. Traces of these reagents may remain on the tissue allografts.

Although this tissue has been tested and screened for relevant communicable diseases and disease agents, human derived tissue allografts may still transmit

infectious agents, known or unknown. Application and use of any allograft tissue may potentially have negative outcomes. Occurrence of complications at the affected site may transpire post-treatment, without early warning signs. These include, but are not limited to: 1) transmission of communicable diseases; 2) transmission of infectious disease agents; and 3) immune rejection of, and/or allergic reaction to the HCT/P. Any adverse outcomes potentially attributable to the HCT/P must be reported promptly to Allotech, LLC.

Caution should be taken when administering this product to immunocompromised individuals, such as patients suffering from HIV or other highly immunocompromised conditions. Although Allotech, LLC has taken great measures to ensure the safety of our allograft products, current technologies cannot preclude the transmission of certain diseases known or unknown.

The end-user shall be held responsible for determining the appropriate application and usage of this product. It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant.

DONOR ELIGIBILITY & SEROLOGICAL TESTING: The donated human tissue utilized for this implant was recovered from a donor screened for risk factors associated with infectious diseases and medical conditions that rule out donation. The donor’s blood was tested for relevant communicable diseases in a laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) or equivalent and registered with the U.S. Food and Drug Administration (FDA) for donor testing. The donor was found negative or non-reactive for:

- HIV-1 / HIV-2 Antibody

- Hepatitis C Virus Antibody
- Hepatitis B Surface Antigen
- Hepatitis B Core Antibody (Total)
- Treponema pallidum (Syphilis)
- Human T-Cell Lymphotropic Virus I/II Antibody
- HIV-1/HCV/ HBV* NAT-TMA

*For donors received after January 01, 2014.

If additional testing was performed (e.g., West Nile Virus), all available test results were reviewed as part of the donor eligibility determination.

A licensed physician determined that the donor met eligibility requirements. The physician utilized available relevant information which may have included but was not limited to donor risk assessment interview, medical/hospital records, donor physical assessment, infectious disease test results, radiology/pathology reports, death certificate and autopsy report (if performed).

MICROBIAL TESTING AND TERMINAL STERILIZATION: Tissue is subjected to microbiological testing and must be free of specific aerobic / anaerobic microorganisms and fungal contaminants whose presence would preclude tissue from processing or transplantation. Following packaging, tissue is terminally sterilized by e-beam irradiation to a minimum sterility assurance level (SAL) of 10⁻⁶.

POTENTIAL COMPLICATIONS / ADVERSE REACTIONS:

- Inherent uncertainty exists in medical and social histories and laboratory testing which may not detect known or unknown pathogens.
- Immune response to transplanted tissue.
- Transmission of known pathogens including Hepatitis B or C, Human T-Cell

Leukemia / Lymphotropic Virus, Human Immunodeficiency Virus 1 & 2, Syphilis, or bacteria. Transmission or causation of diseases of unknown etiology and characteristics.

BEFORE USE: Examine Allograft Packaging and Do Not Use this allograft if: 1) Any of the packaging elements appear to be missing, tampered with or damaged. 2) The product label is severely damaged, illegible, or missing. 3) The expirations date shown on the package label has passed. If any of these conditions exist or are suspected, this allograft should NOT be used.

GRAFT STORAGE: It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant. Store grafts at 15-30°C (59-86°F) in a clean, dry location. Allograft must be stored with proper precautions to prevent contamination.

DIRECTIONS FOR HANDLING AND PREPARATION: Preparation should be performed using aseptic technique. Carefully open the outer pouch and remove sterile inner pouch holding the dermis allograft. The sterile dermis allograft can then be removed from the sterile inner pouch for usage during the surgical procedure. Once the packaging has been opened, the tissue must either be transplanted or discarded. The tissue is intended for use in one patient, on a single occasion only.

Graft Type	Graft Storage	Graft Preparation
Dermis Allograft	Store grafts at 15-30°C (59-86°F) in a clean, dry location.	Cytoplast™ Microderm rapidly hydrates upon contact with blood or saline solution. If pre-hydration is desired, a hydration time of one minute is recommended.

TISSUE TRACKING: The FDA requires patient records to be properly maintained by storing the allograft LOT NUMBER for purposes of tracking the allograft from the donor to the recipient. Please go to our website, <https://osteogenics.com/fda-tissue-tracking> and submit tissue tracking information.

LABELING SYMBOLS: Symbols may be used on package labeling for easy identification.

 Expiration Date	 Single use only	 Lot Number
 Attention, see instruction for use	 Method of Sterilization Using Irradiation	
15°C  30°C Temperature limit 15-30°C (59-86°F)	 Catalog Number	

Red indicator label indicates exposure to irradiation. The tissue’s expiration date is expressed in the format YYYY-MM-DD.

Allotech, LLC processed, stored, and distributed all tissue in accordance with the standards of the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), Health Canada’s Safety of Human Cells, Tissues and Organs for Transplantation Regulations, and applicable State regulations. Osteogenics Biomedical stored and distributed all tissue in accordance with the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), Health Canada’s Safety of Human Cells, Tissues and Organs for Transplantation Regulations, and applicable State regulations.

Responsible Tissue Banks:
Source Establishment, Determines Donor Eligibility & Processing:
RTI Surgical, Inc.
11621 Research Circle
Alachua, FL 32615 USA
CTO Number: 100053
Tel.: +1-800-624-7238

Source Establishment, Processing, Determines Tissue Meets Release Criteria & Distribution:
Allotech, LLC
4620 71st Street, Bldg 78-79
Lubbock, TX 79424
Health Canada CTO Registration Certificate number: 100211
Tel.: +1-806-796-1923

Import and Distribution in Canada:
Nobel Biocare Canada Inc.
150 Via Renzo Drive, Unit 2
Richmond Hill, ON, L4S 0J7
Canada
Tel.: +1-800-939-9394
info.canada@nobelbiocare.com
Canadian reg. CTO no.: 100229

Distribution:
Nobel Biocare USA, LLC
1889 West Mission Boulevard
Pomona, CA 91766
Tel.: +1-800-322-5001
info.usa@nobelbiocare.com

Distribution:
Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424
Tel.: +1-806-796-1923
osteogenics.com

CONTACT INFORMATION: Please contact Allotech at +1-888.796.1923 or +1-806.796.1923 to promptly report any unanticipated or adverse events, or should you require further information. Retain this information for your records.

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