creos™ allo.gain

Allogenic Bone Particulate Package Insert

DESCRIPTION / USE: In accordance with FDA Article 21 CFR Part 1271 and Health Canada's Safety of Human Cells, Tissues and Organs for Transplantation Regulations, this package contains donated human cells, tissues, and cellular and tissue-based products (HCT/Ps). Human allograft bone may be used in a variety of dental reconstructive, regenerative, or periodontal procedures. Tissue is processed and preserved by a variety of techniques, supplied in a range of sizes for surgical use by licensed clinicians and is intended for homologous use only. The tissue is donated human tissue. All tissue is processed and packaged using aseptic technique, freeze dried, and terminally sterilized.

CONTRAINDICATIONS: The presence of infection at the transplantation site is a contraindication for use of musculoskeletal allograft.

WARNINGS AND PRECAUTIONS: Human tissue has the potential to transmit infectious agents. Donor screening, processing treatments and laboratory testing follow stringent specifications to reduce the risk of infectious agent transmission. Do not use if the expiration date has been exceeded or if there is evidence of defects in package or label integrity. Do not sterilize or re-sterilize. It is the responsibility of the clinician to maintain tissue for transplantation according to recommended storage conditions. Do not use if tissue has not been stored according to the recommended STORAGE instructions. This tissue is intended for use by licensed clinicians.

Allografts containing demineralized ground bone are processed using some or all of the following agents: ascorbic acid, hydrochloric acid, hydrogen peroxide, povidone iodine, and sodium phosphate buffer. Allografts containing mineralized ground bone are processed using some or all of the following agents: ascorbic acid, hydrogen peroxide, isopropyl alcohol, and povidone iodine. Traces of these reagents may remain on the tissue allografts.

Although this tissue has been tested and screened for relevant communicable diseases and disease agents, and processed under aseptic conditions, application and use of any allograft tissue may potentially have negative outcomes. Occurrence of complications at the affected site may transpire post-treatment, without early warning signs. These include, but are not limited to: 1) transmission of communicable diseases; 2) transmission of infectious disease agents; and 3) immune rejection of, and/or allergic reaction to the HCT/P. Any adverse outcomes potentially attributable to the HCT/P must be reported promptly to Allotech, LLC.

Caution should be taken when administering this product to immunocompromised individuals, such as patients suffering from HIV or other highly immunocompromised conditions. Although Allotech, LLC has taken great measures to ensure the safety of our allograft products, current technologies cannot preclude the transmission of certain diseases known or unknown. Therefore, Allotech, LLC can make no claims concerning the biological properties and safety of allograft tissue.

The end-user shall be held responsible for determining the appropriate application and usage of this product.

DONOR ELIGIBILITY & SEROLOGICAL TESTING: The donated human tissue utilized for this implant was recovered from a donor screened for risk factors associated with infectious diseases and medical conditions that rule out donation. The donor's blood was tested for relevant communicable diseases in a laboratory certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) or equivalent and registered with the U.S. Food and Drug Administration (FDA) for donor testing. The donor was found negative or non-reactive for:

- HIV-1 / HIV-2 Antibody
- Hepatitis C Virus Antibody
- Hepatitis B Surface Antigen
- Hepatitis B Core Antibody (Total)
- Treponema pallidum (Syphilis)
- Human T-Cell Lymphotropic Virus I/II Antibody
- HIV-1/HCV/ HBV* NAT-TMA

*For donors received after January 01, 2014.

If additional testing was performed (e.g., West Nile Virus), all available test results were reviewed as part of the donor eligibility determination. A licensed physician determined that the donor met eligibility requirements. The physician utilized available relevant information which may have included but was not limited to donor risk assessment interview, medical/hospital records, donor physical assessment, infectious disease test results, radiology/pathology reports, death certificate and autopsy report (if performed).

MICROBIAL TESTING AND TERMINAL STERILIZATION: Prior to packaging, tissue is subjected to microbiological testing at recovery and in the course of processing and must be free of specific aerobic / anaerobic microorganisms and fungal contaminants whose presence would preclude tissue from processing or transplantation. Following packaging, tissue is terminally sterilized by e-beam irradiation to a minimum sterility assurance level (SAL) of 10^{-6}.

POTENTIAL COMPLICATIONS / ADVERSE REACTIONS: Inherent uncertainty exists in medical and social histories and laboratory testing which may not detect known or unknown pathogens. Therefore, the following complications may occur with tissue transplantation: Loss of function or integrity of transplanted tissue with resorption, fragmentation, disintegration, and associated loss of continuity, displacement, bending or fracture. Immune response to transplanted tissue. Transmission of known pathogens including Hepatitis B or C, Human T-Cell Leukemia / Lymphotropic Virus, Human Immunodeficiency Virus 1 & 2, Syphilis, or bacteria. Transmission or causation of diseases of unknown etiology and characteristics.

BEFORE USE: Examine Allograft Packaging. Do Not Use This Allograft If: Any of the packaging elements appear to be missing, tampered with or damaged. The product label is severely damaged, illegible, or missing. The expirations date shown on the package label has passed.

GRAFT STORAGE: It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant. Store freeze dried grafts at 15-30°C (59-86°F) in a clean, dry location. Allograft must be stored with proper precautions to prevent contamination.

DIRECTIONS FOR HANDLING AND PREPARATION: All preparation should be performed using aseptic technique. Carefully open the outer tray of the single blister and aseptically remove the sterile inner container holding the ground bone matrix. The sterile ground bone matrix can then be removed from the sterile inner container for usage during the surgical procedure. Once the packaging has been opened, the tissue must either be transplanted or discarded. The tissue is intended for use in one patient, on a single occasion only.

Graft Type	Graft Storage	Graft Preparation
Ground Bone	Store freeze dried grafts at 15-30°C (59-86°F) in a clean, dry location.	To reconstitute, place graft in sterile basin and cover with sterile isotonic solution for 15-30 minutes, depending on the size of the graft

TISSUE TRACKING: The FDA requires patient records to be properly maintained by storing the allograft LOT NUMBER for purposes of tracking the allograft from the donor to the recipient. Please go to our website and submit tissue tracking information. For United States and in English, visit our website at https://store.nobelbiocare.com/us/en/allograft. For Canada and in English, visit our website at https://store.nobelbiocare.com/ca/en/allograft. For Canada and in French, visit our website at https://store.nobelbiocare.com/ca/fr/allograft.

LABELING SYMBOLS: Symbols may be used on package labeling for easy identification.

Expiration Date	Single use only	LOT Lot Number
Attention, see instruction for use	STERILE R Method of Sterilization Using Irradiation	
15°C 30°C Temperature limit 15-30°C (59-86°F)	REF Catalog Number	

Red indicator label indicates exposure to irradiation. The tissue's expiration date is expressed in the format YYYY-MM-DD.

Allotech, LLC processed, stored, and distributed all tissue in accordance with the standards of the Association for Advancing Tissue and Biologics (AATB), the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), Health Canada's Safety of Human Cells, Tissues and Organs for Transplantation Regulations, and applicable State regulations. Nobel Biocare stored and distributed all tissue in accordance with the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), Health Canada's Safety of Human Cells, Tissues and Organs for Transplantation Regulations, and applicable State regulations.

Responsible Tissue Banks:
Source Establishment, Determines Donor Eligibility & Processing:
RTI Surgical, Inc.
11621 Research Circle
Alachua, FL 32615 USA
Health Canada CTO Registration Certificate Number: 100053
Tel.: +1-800-624-7238

Source Establishment, Processing, Determines Tissue Meets Release Criteria & Distribution:
Allotech, LLC
4620 71st Street, Bldg 78-79
Lubbock, TX 79424, USA
Health Canada CTO Registration Certificate number: 100211
Tel.: 1-806-796-1923

Import and Distribution in Canada:
Nobel Biocare Canada Inc.
150 Via Renzo Drive, Unit 2 Richmond Hill, ON, L4S 0J7
Canada
Tel.: +1-800-939-9394
info.canada@nobelbiocare.com
Canadian reg. CTO no.: 100229

Distribution:
Nobel Biocare USA, LLC
22715 Savi Ranch Parkway, Yorba Linda, CA 92887, USA
Tel.: 1-800-322-5001
Info.usa@nobelbiocare.com

Distribution:
Osteogenics Biomedical, Inc.
5609 58th Street
Lubbock, TX 79424, USA
Tel.: +1-806-796-1923
osteogenics.com

CONTACT INFORMATION: Please contact Allotech at +1-888.796.1923 or +1-806.796.1923 to promptly report any unanticipated or adverse events, or should you require further information. Retain this information for your records.

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