

About Resorba[®]

The Resorba[®] company was founded in September 1931. While new to the US market, Resorba is one of the top 3 suture manufactures in Europe. Resorba[®] dental sutures are made from unique materials that are designed especially for dentistry.

Osteogenics Biomedical is proud to be the exclusive North American distributor of the premium German suture brand.



Glycolon[™] (PGA/PCL)

- Absorbable monofilament
- Maintains 50% tensile strength for 11-13 days
- Sizes range from USP 4/0 to USP 6/0



PGA Resorba[™] (PGA)

- Absorbable multifilament (coated)
- Maintains 50% tensile strength for up to 21 days
- Sizes range from USP 4/0 to USP 6/0




Resolon[™] (Nylon)

- Non-absorbable monofilament
- Heat treated to increase softness and flexibility
- Sizes range from USP 4/0 to USP 7/0



Resolon Twist[™] (Nylon)

- Non-absorbable pseudo-monofilament
- Coated pseudo-monofilament reduces drag
- Sizes range from USP 3/0 to USP 4/0



Available with Black Needles

- ▶ Manufactured from premium tempered black 300 series stainless steel
- ▶ Black steel provides more contrast in the mouth for easier visualization
- ▶ Dark color doesn't reflect light, which is advantageous when suturing under magnification

